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Daily checklist

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- Quick release locked ✓
- Head rest is in correct height ✓
- Belt is tightened and adjusted ✓
- Foot rest is in correct height ✓
- Foot straps are tightened ✓
- Check frequently if bolts and screws are tightened ✓

goguard by Bikevision ApS
+45 8282 3232
info@goguard.dk
www.goguard.eu

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*The chair comes with a separate seat cover of your choice

Mount the cover of your choice by pulling the safety belt through the holes. The cover is fastened by five push buttons.

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Mount the neck support in the two holes at the top of the chair. Lock the neck support in the chosen position by pressing the 2 levers together on the backside of the chair.

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Assemble the safety belt by “flexing” the two belt parts and then press the belt clips in the lock. To open, remember to press both sides while pushing the big central button. Remember to adjust the length of the belt to the child’s size.

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Leg supports l+r are pressed together around the center tube and fastened with the two following screws. Keep pressing the parts towards each other until the holes are visible.

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Leg supports can be adjusted by pressing the two hooks slightly together . Adjust as desired.

The 2 parts of the rear bracket are fastened with 4 bolts around the rear stay.

- The front bracket is fastened with 2 bolts around the saddle tube.

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Push the monotube downwards through the rear bracket with the quick release being open. Push in the safety button. The bracket should be adjusted in the correct height, so the monotube is fit for a possible change of the back position. Bracket is fastened in position with a key. The center of the chair should follow the centre of the bike. When in position, close the quick release and tighten. The chair should not be able to move.

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The rear support stays can be fastened by pulling them all the way back. The bracket should be positioned so the support stays are as vertical as possible. When in position, tighten.

Lift the chair horizontally, 5 cm above the luggage carrier. Then tighten the bolt with an 18mm tool.

The chair is now ready. Remember you can adjust to rest position by opening the quick release and move the monotube. Always remember to lock again.